

## SCIENTIFIC BOOKLET





## 10 scientific references

1.	Gait-like vibration training improves gait abilities: a case report of a 62-year-old person with a chronic incomplete spinal cord injury - Barthelemy, 2016	€
2.	Alternate rhythmic vibratory stimulation of trunk muscles affects walking cadence and velocity in Parkinson's disease - De Nunzio, 2009	€
З.	Vibration Elicits Involuntary, Step-Like Behavior in Individuals With Spinal Cord Injury - Field-Fote, 2012	€
4.	Effect of illusory kinesthesia on hand function in patients with distal radius fractures: a quasi-randomized controlled study - Imai, 2017	€
5.	Effect of a local vibration stimulus training programme on postural sway and gait in chronic stroke patients - Lee, 2013	€
6.	Focal vibration in neurorehabilitation - Murillo, 2014	€
7.	Illusory movements prevent cortical disruption caused by immobilization - Roll, 2012	€
8.	Segmental muscle vibration improves reaching movement in patients with chronic stroke. A randomized controlled trial - Tavernese, 2013	€
9.	Short-Term Effects of Focal Muscle Vibration on Motor Recovery After Acute Stroke - Toscano, 2019	€
10.	Effectiveness of matrix-rhythm therapy on increased muscle tone, balance and gait parameters in stroke survivors - Unal, 2020	€



#### **Mode 1 | Mobility** Functional Proprioceptive Stimulations

#### **Clinical Benefits**

- FPS applied very early during the acute phase before passive mobilization or conventional therapy
- + FPS preserve sensory-motor interaction
- + FPS can compensate for the lack of feedback due to long-term immobility
- + FPS activate the cortical networks of active movements guiding cortical plasticity
- + FPS reduce the side-effects of immobility and the need for re-training
- + FPS trigger motor responses consistent with the sensation evoked
- + FPS enhance the coordination and limit the co-contractions
- + FPS accelerate the recovery of mobility and motor control
- + FPS reduce the duration of the treatment while improving patient comfort
- + Positive additive effect (enhance the corticomotor excitability) of the combined work of FPS/motor imagery/passive manipulation or active participation





#### Mode 2 | Antalgy

## 10 scientific references

1.	After-effects of neck muscle vibration on sensorimotor function and pain in neck pain patients and healthy controls - Beinert, 2019	•
2.	The analgesic effect of localized vibration: a systematic review - Casale, 2022	€
З.	Proprioceptive feedback enhancement induced by vibratory stimulation in complex regional pain syndrome type I - Gay, 2007	€
4.	Substance P-like immunoreactivity and analgesic effects of vibratory stimulation on patients suffering from chronic pain - Guieu, 1993	€
5.	How does vibration reduce pain? - Hollins, 2014	€
6.	Mechanisms of pain relief by vibration and movement - Kakigi, 1992	•
7.	Pain alleviation by vibratory stimulation - Lundeberg, 1984	€
8.	Somatosensory rehabilitation for allodynia in complex regional pain syndrome of the upper limb - Packham, 2018	€
9.	Static mechanical allodynia (SMA) is a paradoxical painful hypo-aesthesia - Spicher, 2008	€

10. Vibration reduces thermal pain in adjacent dermatomes - Yarnitsky, 1997



 $\overline{\mathbf{G}}$ 

#### Mode 2 | Antalgy

#### **Clinical benefits**

- + Reduction of pain through the stimulation of skin mechanoreceptors, at distance to the painful area, through the Gate Control effect: pain can be reduced by simultaneous activation of nerve fibres that conduct non-noxious stimuli
- + Make the patient gradually bear increasingly significant dose of mechanical vibrations
- + Application of FV with light pressure generates a powerful sensory flow transmitted to the nerve centers via large-diameter afferent fibers to give rise painless tactile sensations
- + Low-intensity mechanical stimuli, not aggressive enough to stimulate nociceptors, activates the inhibiting interneurons: closed-door theory of Melzack & Wall
- + Reduction of pain and improvement of comfort allow to do more coventional therapy
- + Inhibitory effects of peripheral origin would combine with those of central origin, resulting from spinal or cortical treatment of painful informations
- Possible desensitization : proximal to distal examination with focal vibrations in zig zag along the path of the nerve
- + Allodynia : remote vibrotactile counter-stimulation at distance to the painful area
- + Preserving one's maximum force production capacities allow to carry out daily activities but also to prevent osteoarthritis
- + In case of mechanical allodynia, it is crucial to use FV to ensure that the correct diagnosis is made before starting treatment





## 10 scientific references

1.	Localized muscle vibration reverses quadriceps muscle hypotrophy - Benedetti, 2017	€
2.	Increases in muscle activity produced by vibration of the thigh muscles during locomotion in chronic human SCI - Cotey, 2009	€
3.	Is the Focal Muscle Vibration an Effective Motor Conditioning Intervention ? - Fattorini, 2021	€
4.	Potentiation of muscle strength by focal vibratory stimulation on quadriceps femoris - Feltroni, 2018	€
5.	Improvement of Stance Control and Muscle Performance Induced by Focal Muscle Vibration - Filipi, 2009	€
6.	Effect of Segment-Body Vibration on Strength Parameters - Goebel, 2015	€
7.	Effect of vibratory stimulation training on maximal force and flexibility - Issurin, 2008	€
8.	Effects of focal muscle vibration on physical functioning in patients with knee osteoarthritis: a randomized controlled trial - Rabini, 2015	€
9.	Facilitation of triceps brachii muscle contraction by tendon vibration after chronic cervical spinal cord injury - Ribot-Ciscar, 2003	€
10.	Influence of vibration on endurance of maximal isometric contraction - Samuelson, 1989	€



## **Clinical benefits**

- + Long-lasting positive motor conditioning effect
- + Increase muscle tone
- + Limit muscle hypotrophy due to immobility or lack of movement
- + Strenghten motor function
- + Carry out important proprioceptive work
- + Additional positive effect of isometric contractions of the stimulated muscles
- + Large acceptance, confortable sensations and no negative collateral effect
- + Increase balance: plasticity of the motor system can be enhanced by proprioceptive inputs
- + Increasing muscle activity allow to decrease the amount of required support/ assistance & improve the efficacy of the treatment
- + Optimal effects with multiple sessions of treatment (4-5 sessions/week or more)





#### Mode 4 | Spasticity

## 10 scientific references

1.	Effectiveness of Focal Muscle Vibration on Hemiplegic Upper Extremity Spasticity in Individuals With Stroke: A Systematic Review - Alashram, 2019	<b>Э</b>
2.	Focal muscle vibration in the treatment of upper limb spasticity - Caliandro, 2012	€
З.	Localized 100 Hz vibration improves function and reduces upper limb spasticity: a double- blind controlled study - Casale, 2014	€
4.	Preliminary evidence of focal muscle vibration effects on spasticity due to cerebral palsy in a small sample of Italian children - Celetti, 2011	€
5.	Brain oscillatory activity correlates with the relief of post-stroke spasticity following focal vibration - Li, 2022	€
6.	Decrease of spasticity with muscle vibration in patients with spinal cord injury - Murillo, 2011	€
7.	Anti-spastic effects of the direct application of stimuli to the spastic muscles - Noma, 2012	€
8.	Vibration therapy in patients with cerebral palsy a systematic review - Ritzmann, 2018	€
9.	The effects of robotic gait neurorehabilitation and focal vibration combined treatment in adult cerebral palsy - Rutovic, 2019	€
10.	Motor Recovery After Stroke: From a Vespa Scooter Ride Over the Roman Sampietrini to Focal Muscle Vibration (fMV) Treatment - Toscano, 2020	€



#### Mode 4 | Spasticity

#### **Clinical benefits**

- + Safe & well tolerated ; easy to perform at beside from the acute phase to the chronic one
- + Long-lasting regulation of the muscle tone between agonist & antagonist
- + Can preserve muscular architecture
- + Focal vibration combined with physiotherapy is better and faster than physiotherapy alone in controlling spasticity and improving motor function
- + Different type of application depending on the site and type/level of spasticity
- + Repeated muscle vibration produces a repeated sensory input that reaches, via la fiber affetent input, S1/M1 cortical areas and induce plasticity
- + FV can provoke cortical and spinal plasticity leading to the reduction of spasticity
- + Increase excitability in the primary motor cortex especially in S1-M1: biomarker of intrinsic plasticity-related mechanism for the reduction of spasticity
- + The reduction of hypertonia leads to less pain and allow the patients to do more conventional therapy or robotic one
- + Possible reduction of the frequency and amount of medication or toxin taken by the patient









# SENSORY-MOTOR THERAPY



Vibramoov Physio

Vibramoov Physio is a Medical Device, CE marked (Medical Device European Regulation) designed and manufactured upon ISO 13485: 2016 standard



25 Place de l'Encas | ZI Saint-Maurice | 04100 Manosque | France +33(0)4 92 79 08 56 www.technoconcept.com